



Laerdal
helping save lives



Laerdal Silicone Resuscitator

Reusable BVM

The reusable Laerdal Silicone Resuscitator (LSR) enables medical personnel to facilitate assisted ventilation of patients in pre-hospital, in-hospital and post-hospital care. The LSR is designed for convenience, long product life and meets the needs of both the rescuer and the patient. The LSR comes in adult, child and preterm sizes.



Available in carton, compact case and display case



Full range of mask sizes



Available in adult, child and pre-term sizes

Laerdal Silicone Resuscitator

Reusable BVM

KEY FEATURES

- Easy assembly
- Standard 15/22mm connectors
 - HMEF filters with a 15M/22F connector will fit on the LSR
- Designed to prevent rebreathing of exhaled air
- 35cm H₂O Pressure Relief Valve reduces risk of excessive ventilation pressure
- Long product life
- Can be high-level disinfected and sterilized
- Efficient integrated Intake/Reservoir Valve with reservoir bag ensures oxygen economy and increased oxygen concentration
- Foldable for compact storage to save space
- Available in three packaging configurations: carton, compact case and display case
- Easy decontamination minimizes cleaning time and reduces inventory needs
- Can be reused for years (with proper care)

MODELS, CONTENT AND PACKAGING



	ADULT (Patients >55 lbs.)			PEDIATRIC (Patients 5.5 to 55lbs.)				PRETERM (Patients < 5.5 lbs.)	
	Basic w/o Mask	Standard	Complete	Basic w/o Mask	Standard Child	Standard Term-Baby	Complete	Basic w/o Mask	Complete
Carton	870050	870052	870051	860050	860052	860056	860051	850050	850051
Compact Case			870053				860053		850053
Display Case			870055				860055		850055
Bag 1600 ml	•	•	•						
500 ml				•	•	•	•		
240 ml								•	•
Pt Valve w/o Relief Valve	•	•	•						
w/ Relief Valve				•	•	•	•	•	•
Mask 0/0									•
0/1						•	•		•
2							•		
3/4 w/Mask Cover			•		•		•		
4/5+ w/Mask Cover		•	•						
Intake/Reservoir Valve	•	•	•	•	•	•	•	•	•
Reservoir Bag 2600	•	•	•						
600				•	•	•	•	•	•
DFU Instructions for Use	•	•	•	•	•	•	•	•	•



Available in carton, compact case and display case



Full range of mask sizes



Available in adult, child and pre-term sizes

TECHNICAL SPECIFICATIONS

Temperature and Environment	
Operating Temperature Humidity	-18°C to 60°C (0°F to 140°F) 15% to 95% rH
Long term Storage Temperature Humidity	40°C to 70°C (-40°F to 158°F) 40% to 95% rH
Adult Model	
Ventilation Bag volume	1600 ml
Reservoir Bag volume	2600 ml
Child Model	
Ventilation Bag volume	500 ml
Reservoir Bag volume	600 ml
Preterm Model	
Ventilation Bag volume	240 ml
Reservoir Bag volume	600 ml
Attainable delivery volume	
Adult	Approximately 800 ml
Child	Approximately 320 ml
Preterm	Approximately 150 ml



MASKS

Catalogue #	Description
851500	LSR Silicone mask No.00
851600	LSR Silicone mask No.0/1
851700	LSR Silicone mask No.2
860220	Child Silicone mask 3-4 w/Multi Function Mask Cover
870220	Adult Silicone mask 4-5+ w/Multi Function Mask Cover
860221	Child Silicone mask 3-4 w/o Multi Function Mask Cover
870221	Adult Silicone mask 4-5+ w/o Multi Function Mask Cover
872220	Adult & Child silicone Mask w/Multi Function Mask Covers

SPARE PARTS/CONSUMABLES

Catalogue #	Description
540103	LSR Lip valve
851252	Pressure Relief Valve 35 cm H2O
851103	LSR Lock clips,(pkg-10)
540105	LSR Disk membranes, pkg. 10
511700	LSR Wall bracket
531904	LSR Reservoir valve lid
560200	LSR Patient Valve
871300	Attachment ring for head strap
850150	Preterm Bag, 240 ml
860150	Pediatric Bag, 500 ml
870150	Adult Bag, 1600 ml
860410	LSR Compact Case
875400	Intake Reservoir Valve
531901	LSR O ² reservoir 2.6 litre
531906	LSR O ² reservoir 2.6 liters x 50
551901	LSR O ² reservoir 0.6 litre
510404	LSR Intake membranes pkg. 10
871950	Flap Valves 2 pk.
510103	Cap for LSR Intake Valve spare part pkg. 3
860420	LSR Compact Case Adult
870600	LSR Display case cpl. Adult
860300	Display Case, Pediatric
850700	Display Case, Preterm
530400	Airways, set of 4
871950	Flap Valves 2 pk.
870708	Disposable O ² -reservoir bag, 2600ml, pkg.20
851250	Patient Valve w/35cm H ₂ O Pressure Relief Valve
521100	Wall Mount Pediatric/ Preterm Display Case

ACCESSORIES

Catalogue #	Description
871000	Silicone extension tube, 28 cm
850900	LSR Manometer connector
850500	Expiration diverter (OD 30 mm)
531907	LSR Intake Valve Outer Part (23 mm OD)
865200	Multi Function Mask Cover for Mask 3-4
875200	Multi-Function Mask Cover 4-5+
870400	LSR Head strap
870120	LSR Hanging loop
572000	Wall mount, Adult Display Case



Laerdal
helping save lives